


MRGPRX1 agonist 2

Chemical Properties

CAS No. :	302952-43-8
Formula:	C ₁₅ H ₁₄ N ₂ O ₂ S
Molecular Weight:	270.35
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	MRGPRX1 agonist 2 (compound 1a) is a potent positive allosteric modulator of Mas-related G protein-coupled receptor X1 (MRGPRX1) with an EC ₅₀ value of 0.48 μM, and has potential for use in neuropathic pain research [1].
Targets(IC50)	Others,MRGPR

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6989 mL	18.4945 mL	36.9891 mL
5 mM	0.7398 mL	3.6989 mL	7.3978 mL
10 mM	0.3699 mL	1.8495 mL	3.6989 mL
50 mM	0.074 mL	0.3699 mL	0.7398 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481