

5-[5-(methoxycarbonyl)-5,8a-dimethyl-2-methylidene-decahydronaphthalen-1-yl]-3-methylpentanoic acid

Chemical Properties

CAS No. :

Formula: C₂₁H₃₄O₄

Molecular Weight: 350.5

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	5-[5-(methoxycarbonyl)-5,8a-dimethyl-2-methylidene-decahydronaphthalen-1-yl]-3-methylpentanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8531 mL	14.2653 mL	28.5307 mL
5 mM	0.5706 mL	2.8531 mL	5.7061 mL
10 mM	0.2853 mL	1.4265 mL	2.8531 mL
50 mM	0.0571 mL	0.2853 mL	0.5706 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481