

Histone H2BK12ac (1-22)-GGK-biotin amide (trifluoroacetate salt)

**Chemical Properties**

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

**Biological Description**

Description	Histone H2BK12ac (1-22)-GGK-biotin amide (trifluoroacetate salt) can be used in related research in the field of life sciences, and its product number is T36095.
Targets(IC50)	Others

**Solubility Information**

Solubility	DMSO: Soluble (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481