

Compound N094-0016

Chemical Properties


CAS No. :

Formula: C₁₅H₁₇NO₄

Molecular Weight: 275.3

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

| | |
|---------------|--|
| Description | Compound N094-0016 is a useful organic compound for research related to life sciences and the catalog number is T131160. |
| Targets(IC50) | Others |

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|-----------|-----------|
| 1 mM | 3.6324 mL | 18.162 mL | 36.324 mL |
| 5 mM | 0.7265 mL | 3.6324 mL | 7.2648 mL |
| 10 mM | 0.3632 mL | 1.8162 mL | 3.6324 mL |
| 50 mM | 0.0726 mL | 0.3632 mL | 0.7265 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481