


## Praziquantel-D11

## Chemical Properties

CAS No. :	1246343-36-1
Formula:	C <sub>19</sub> H <sub>24</sub> N <sub>2</sub> O <sub>2</sub>
Molecular Weight:	323.47
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Praziquantel-D11 is the deuterium labeled Praziquantel (TMSM-1935), and is an agent of anthelmintic.
Targets(IC50)	Parasite

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0915 mL	15.4574 mL	30.9148 mL
5 mM	0.6183 mL	3.0915 mL	6.183 mL
10 mM	0.3091 mL	1.5457 mL	3.0915 mL
50 mM	0.0618 mL	0.3091 mL	0.6183 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481