


## Ammonium, (4-(p-trimethylammoniumphenyl)butyl)trimethyl-, diiodide

## Chemical Properties

CAS No. :	63951-20-2
Formula:	C <sub>16</sub> H <sub>30</sub> I <sub>2</sub> N <sub>2</sub>
Molecular Weight:	504.23
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Ammonium, (4-(p-trimethylammoniumphenyl)butyl)trimethyl-, diiodide is a bioactive chemical.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9832 mL	9.9161 mL	19.8322 mL
5 mM	0.3966 mL	1.9832 mL	3.9664 mL
10 mM	0.1983 mL	0.9916 mL	1.9832 mL
50 mM	0.0397 mL	0.1983 mL	0.3966 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481