

Isotretinoin

Chemical Properties

CAS No. : 4759-48-2


Formula: C₂₀H₂₈O₂

Molecular Weight: 300.44

Storage: The compound is unstable in solution. Please use soon

Powder: -20°C for 3 years

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Isotretinoin (13-cis-Retinoic acid) binds to and activates nuclear retinoic acid receptors (RARs), which act as transcription factors that promote cell differentiation and apoptosis. This naturally occurring retinoic acid has potential antineoplastic activity, exhibits immunomodulatory and anti-inflammatory responses, and inhibits ornithine decarboxylase, thereby decreasing polyamine synthesis and keratinization.
Targets(IC50)	Retinoid Receptor,Endogenous Metabolite,Autophagy
In vitro	Isotretinoin directly interferes with the development of cranial neural crest cells. [1] Isotretinoin selectively affects neural crest cells by decreasing their cell-substratum adhesion. [2]
In vivo	Isotretinoin (500 ng/mL) and its main metabolite in the human, 4-oxo-isotretinoin, induce malformations similar to those seen in vivo. [1] Isotretinoin impairs explicit memory in Stage 2, but retention tests one month after Isotretinoin exposure ended, indicated recovery from this explicit memory impairment and evidence of enhanced implicit memory in the 10 mg and 15 mg ISO rats. [3] Isotretinoin slows the recovery of rod signaling after exposure to an intense bleaching light, and that rhodopsin regeneration is markedly slowed. Isotretinoin is also found to protect rat photoreceptors from light-induced damage. [4] Isotretinoin blocks the formation of A2E biochemically and the accumulation of lipofuscin pigments by electron microscopy. Isotretinoin also blocks the slower, age-dependent accumulation of lipofuscin in wild-type mice. [5]

Solubility Information

Solubility	Ethanol: 1 mg/mL (3.33 mM),Heating is recommended. (The compound is unstable in solution, please use soon.) DMSO: 250 mg/mL (832.11 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 2 mg/mL (6.66 mM),Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may</i>

A DRUG SCREENING EXPERT

In vivo Formulation	<i>vary and should be modified based on specific experimental conditions.</i>
---------------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3285 mL	16.6423 mL	33.2845 mL
5 mM	0.6657 mL	3.3285 mL	6.6569 mL
10 mM	0.3328 mL	1.6642 mL	3.3285 mL
50 mM	0.0666 mL	0.3328 mL	0.6657 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Webster WS, et al. J Craniofac Genet Dev Biol, 1986, 6(3), 211-222.

Li T, Lu D, Yao C, et al. Kansl1 haploinsufficiency impairs autophagosome-lysosome fusion and links autophagic dysfunction with Koolen-de Vries syndrome in mice. Nature Communications. 2022, 13(1): 1-16.

Qiu Y, Sun Y, Xu D, et al. Screening of FDA-approved drugs identifies sutent as a modulator of UCP1 expression in brown adipose tissue. EBioMedicine. 2018, 37: 344-355.

Smith-Thomas L, et al. Dev Biol, 1987, 123(1), 276-281.

Dopheide MM, et al. Pharmacol Biochem Behav, 2008, 91(2), 243-251.

Sieving PA, et al. Proc Natl Acad Sci U S A, 2001, 98(4), 1835-1840.

Radu RA, et al. Proc Natl Acad Sci U S A, 2003, 100(8), 4742-4747.

Qiu Y, Sun Y, Xu D, et al. Screening of FDA-approved drugs identifies sutent as a modulator of UCP1 expression in brown adipose tissue[J]. EBioMedicine. 2018, 37: 344-355.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481