


5-Methoxy-1,7-diphenyl-3-heptanone

Chemical Properties

CAS No. : 100667-54-7
Formula: C₂₀H₂₄O₂
Molecular Weight: 296.4
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

| | |
|-------------|--|
| Description | 5-Methoxy-1,7-diphenyl-3-heptanone is a natural product. |
|-------------|--|

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|------------|------------|
| 1 mM | 3.3738 mL | 16.8691 mL | 33.7382 mL |
| 5 mM | 0.6748 mL | 3.3738 mL | 6.7476 mL |
| 10 mM | 0.3374 mL | 1.6869 mL | 3.3738 mL |
| 50 mM | 0.0675 mL | 0.3374 mL | 0.6748 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481