

Compound T129990

Chemical Properties


CAS No. :

Formula: C36H56O14

Molecular Weight: 712.83

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	4-[[[(1R,3aS,7R,9R,9aS,11aR)-7-[[[(2R)-5-[[[3,5-dihydroxy-4-methoxy-6-methyloxan-2-yl]oxy]-3,4-dihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy]-3a,9-dihydroxy-9a,11a-dimethylhexadecahydro-1H-cyclopenta[a]phenanthren-1-yl]-2,5-dihydrofuran-2-one]] is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4029 mL	7.0143 mL	14.0286 mL
5 mM	0.2806 mL	1.4029 mL	2.8057 mL
10 mM	0.1403 mL	0.7014 mL	1.4029 mL
50 mM	0.0281 mL	0.1403 mL	0.2806 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481