

(2E)-4-hydroxydodec-2-enedioic acid

Chemical Properties

CAS No. :

Formula: C₁₂H₂₀O₅

Molecular Weight: 244.29

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	(2E)-4-hydroxydodec-2-enedioic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.0935 mL	20.4675 mL	40.935 mL
5 mM	0.8187 mL	4.0935 mL	8.187 mL
10 mM	0.4093 mL	2.0467 mL	4.0935 mL
50 mM	0.0819 mL	0.4093 mL	0.8187 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481