


## N-methyl-(4-phenyltetrahydropyran-4-yl)methylamine

## Chemical Properties

CAS No. :	958443-30-6
Formula:	C <sub>13</sub> H <sub>19</sub> NO
Molecular Weight:	205.3
Storage:	Store under nitrogen Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	N-methyl-(4-phenyltetrahydropyran-4-yl)methylamine, with CAS No. 958443-30-6, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. N-methyl-(4-phenyltetrahydropyran-4-yl)methylamine provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.8709 mL	24.3546 mL	48.7092 mL
5 mM	0.9742 mL	4.8709 mL	9.7418 mL
10 mM	0.4871 mL	2.4355 mL	4.8709 mL
50 mM	0.0974 mL	0.4871 mL	0.9742 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481