


## Tyrphostin AG17

## Chemical Properties

CAS No. :	71308-35-5
Formula:	C <sub>12</sub> H <sub>10</sub> N <sub>2</sub> O
Molecular Weight:	198.22
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Tyrphostin AG17 is a biochemical that promotes adipocyte apoptosis in vivo and is an effective modulator of adipocyte differentiation and WAT hypertrophy in vitro and in vivo. It is thought to possibly be useful as a pharmacological obesity treatment. It is a tyrosine kinase antagonist that has also been found to inhibit tumor cell growth.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.0449 mL	25.2245 mL	50.449 mL
5 mM	1.009 mL	5.0449 mL	10.0898 mL
10 mM	0.5045 mL	2.5224 mL	5.0449 mL
50 mM	0.1009 mL	0.5045 mL	1.009 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481