


JCC76

Chemical Properties

CAS No. : 1021926-22-6
 Formula: C₂₄H₃₂N₂O₄S
 Molecular Weight: 444.59
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	JCC76 is a selective inhibitor of heat shock protein 27 (Hsp27). JCC76 inhibits the proliferation of human epidermal growth factor 2 (Her2).
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2493 mL	11.2463 mL	22.4926 mL
5 mM	0.4499 mL	2.2493 mL	4.4985 mL
10 mM	0.2249 mL	1.1246 mL	2.2493 mL
50 mM	0.045 mL	0.2249 mL	0.4499 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481