

CBN209350

Chemical Properties


CAS No. : 726201-10-1

Formula: C₁₃H₁₂N₄O₃

Molecular Weight: 272.26

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	CBN209350 is a KDM4 inhibitor.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.673 mL	18.3648 mL	36.7296 mL
5 mM	0.7346 mL	3.673 mL	7.3459 mL
10 mM	0.3673 mL	1.8365 mL	3.673 mL
50 mM	0.0735 mL	0.3673 mL	0.7346 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481