


[Ala11,22,28]VIP

Chemical Properties

CAS No. : 291524-04-4
Formula: C139H231N43O39S
Molecular Weight: 3160.7
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	[Ala11,22,28]VIP can be used for related research in the field of life sciences. Its product number is T36636 and CAS number is 291524-04-4.
Targets(IC50)	Others

Solubility Information

Solubility	H2O: Soluble (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.3164 mL	1.5819 mL	3.1639 mL
5 mM	0.0633 mL	0.3164 mL	0.6328 mL
10 mM	0.0316 mL	0.1582 mL	0.3164 mL
50 mM	0.0063 mL	0.0316 mL	0.0633 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481