


PPAR γ agonist 3

Chemical Properties

CAS No. : 2011801-48-0
 Formula: C₂₄H₂₃N₃O
 Molecular Weight: 369.46
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	PPAR γ agonist 3, also known as Compound 18a, is a potent and selective agonist of PPAR γ . This compound does not exhibit cytotoxicity towards both non-resistant and resistant cells. Notably, PPAR γ agonist 3 demonstrates antitumor efficacy exclusively when co-administered with Imatinib [1].
Targets(IC50)	Others,PPAR

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7067 mL	13.5333 mL	27.0665 mL
5 mM	0.5413 mL	2.7067 mL	5.4133 mL
10 mM	0.2707 mL	1.3533 mL	2.7067 mL
50 mM	0.0541 mL	0.2707 mL	0.5413 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481