

Compound N064-0037

Chemical Properties


CAS No. :

Formula: C₃₁H₄₆O₄

Molecular Weight: 482.705

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound N064-0037 is a useful organic compound for research related to life sciences and the catalog number is T131205.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0716 mL	10.3582 mL	20.7164 mL
5 mM	0.4143 mL	2.0716 mL	4.1433 mL
10 mM	0.2072 mL	1.0358 mL	2.0716 mL
50 mM	0.0414 mL	0.2072 mL	0.4143 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481