

ZEN-2759

Chemical Properties

CAS No. : 1616400-50-0


Formula: C17H16N2O2

Molecular Weight: 280.32

Store at low temperature

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	ZEN-2759 is a potent BET small molecule inhibitor of BRD4 (BD1), BRD4 (BD2), and BRD4 (BD1BD2).
Targets(IC50)	Epigenetic Reader Domain

Solubility Information

Solubility	DMSO: 60 mg/mL (214.04 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 2 mg/mL (7.13 mM), Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.5674 mL	17.8368 mL	35.6735 mL
5 mM	0.7135 mL	3.5674 mL	7.1347 mL
10 mM	0.3567 mL	1.7837 mL	3.5674 mL
50 mM	0.0713 mL	0.3567 mL	0.7135 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Kharenko OA, et al. Design and Characterization of Novel Covalent Bromodomain and Extra-Terminal Domain (BET) Inhibitors Targeting a Methionine. J Med Chem. 2018 Sep 27;61(18):8202-8211.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481