

SCH-79687

Chemical Properties


CAS No. : 224585-45-9

Formula: C₁₈H₁₆Cl₂N₄O

Molecular Weight: 375.25

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	SCH-79687 is a histamine antagonist selective for the H3 subtype. This H3 antagonist may play a role in allergic rhinitis.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6649 mL	13.3245 mL	26.6489 mL
5 mM	0.533 mL	2.6649 mL	5.3298 mL
10 mM	0.2665 mL	1.3324 mL	2.6649 mL
50 mM	0.0533 mL	0.2665 mL	0.533 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481