

m-PEG-acrylate (MW 10000)

**Chemical Properties**

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

**Biological Description**

Description	m-PEG-acrylate (MW 10000) is a polyethylene glycol (PEG)-based linker used in the synthesis of Proteolysis Targeting Chimeras (PROTACs)[1].
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two ligands linked together—one targeting an E3 ubiquitin ligase and the other targeting the desired protein—utilizing the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

## Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481