


Physalin L

Chemical Properties

CAS No. :	113146-74-0
Formula:	C ₂₈ H ₃₂ O ₁₀
Molecular Weight:	528.55
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Physalin L is a natural product, it shows a distinct fluorescence spot under UV 365 nm with good separation.
Targets(IC50)	Endogenous Metabolite, NO Synthase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.892 mL	9.4598 mL	18.9197 mL
5 mM	0.3784 mL	1.892 mL	3.7839 mL
10 mM	0.1892 mL	0.946 mL	1.892 mL
50 mM	0.0378 mL	0.1892 mL	0.3784 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Physalins with anti-inflammatory activity are present in *Physalis alkekengi* var. *franchetii* and can function as Michael reaction acceptors. *Steroids*. 2012 Apr;77(5):441-7.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481