


PROTAC IRAK4 degrader-6

Chemical Properties

CAS No. :	2360530-72-7
Formula:	C42H41F3N12O8
Molecular Weight:	898.861
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	PROTAC IRAK4 degrader-6 is a potent Cereblon-based compound designed to degrade [interleukin-1 receptor-associated kinase 4 (IRAK4)].
Targets(IC50)	Others,IRAK,PROTACs
In vitro	PROTAC IRAK4 Degrader-6 is a PROTAC-based agent designed to target and degrade interleukin-1 receptor-associated kinase 4 (IRAK4)[1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1125 mL	5.5626 mL	11.1252 mL
5 mM	0.2225 mL	1.1125 mL	2.225 mL
10 mM	0.1113 mL	0.5563 mL	1.1125 mL
50 mM	0.0223 mL	0.1113 mL	0.2225 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

https://scifinder-n.cas.org/patent-viewer?docuri=g__2-GF3YbbKfwwYAYLNmYDipK1sFeQGP9cmlPB_4&markedFullTextKey=GFB6kc_cFl60fBy8t-u8qupJE998gaydSFApjgE-p5c.pdf

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481