

Compound 0449-0135

Chemical Properties


CAS No. :

Formula: C₄₁H₆₄O₁₃

Molecular Weight: 764.95

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound 0449-0135 is a useful organic compound for research related to life sciences and the catalog number is T131647.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3073 mL	6.5364 mL	13.0727 mL
5 mM	0.2615 mL	1.3073 mL	2.6145 mL
10 mM	0.1307 mL	0.6536 mL	1.3073 mL
50 mM	0.0261 mL	0.1307 mL	0.2615 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481