


## 21-hydroxy Eplerenone

## Chemical Properties

CAS No. :	334678-67-0
Formula:	C <sub>24</sub> H <sub>30</sub> O <sub>7</sub>
Molecular Weight:	430.497
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	21-hydroxy Eplerenone is a major metabolite of the mineralocorticoid receptor antagonist eplerenone. It is formed from eplerenone by the cytochrome P450 (CYP) isoform CYP3A4.
Targets(IC50)	Glucocorticoid Receptor, Others

## Solubility Information

Solubility	Chloroform: Slightly soluble Methanol: Slightly soluble ( $< 1$ mg/ml refers to the product slightly soluble or insoluble)
------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3229 mL	11.6144 mL	23.2288 mL
5 mM	0.4646 mL	2.3229 mL	4.6458 mL
10 mM	0.2323 mL	1.1614 mL	2.3229 mL
50 mM	0.0465 mL	0.2323 mL	0.4646 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481