


## Dorrigocin B

## Chemical Properties

CAS No. :	158446-30-1
Formula:	C <sub>27</sub> H <sub>41</sub> N <sub>0</sub> O <sub>8</sub>
Molecular Weight:	507.624
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Dorrigocin B is a natural product that can be used as a reference standard. The CAS number of Dorrigocin B is 158446-30-1.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.970 mL	9.8499 mL	19.6998 mL
5 mM	0.394 mL	1.970 mL	3.940 mL
10 mM	0.197 mL	0.985 mL	1.970 mL
50 mM	0.0394 mL	0.197 mL	0.394 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481