


## Bornaprine

## Chemical Properties

CAS No. :	20448-86-6
Formula:	C <sub>21</sub> H <sub>31</sub> NO <sub>2</sub>
Molecular Weight:	329.48
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Bornaprine is an anticholinergic drug used to treat Parkinson's disease due to inhibition of muscarinic or nicotinic transmission.
Targets(IC50)	Others,Cholinesterase (ChE)

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0351 mL	15.1754 mL	30.3509 mL
5 mM	0.607 mL	3.0351 mL	6.0702 mL
10 mM	0.3035 mL	1.5175 mL	3.0351 mL
50 mM	0.0607 mL	0.3035 mL	0.607 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481