


## CDKi Hydrochloride

## Chemical Properties

CAS No. :	70035-83-5
Formula:	C <sub>4</sub> H <sub>6</sub> N <sub>4</sub> O
Molecular Weight:	126.12
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	CDKi hydrochloride is a potent and selective ATP competitive inhibitor of cycB-CDK1, cycE-CDK2, and cycD1-CDK4.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.929 mL	39.6448 mL	79.2896 mL
5 mM	1.5858 mL	7.929 mL	15.8579 mL
10 mM	0.7929 mL	3.9645 mL	7.929 mL
50 mM	0.1586 mL	0.7929 mL	1.5858 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481