

Compound N016-0045

Chemical Properties


CAS No. :

Formula: C43H68N4O14

Molecular Weight: 865.031

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound N016-0045 is a useful organic compound for research related to life sciences and the catalog number is T131339.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.156 mL	5.7801 mL	11.5603 mL
5 mM	0.2312 mL	1.156 mL	2.3121 mL
10 mM	0.1156 mL	0.578 mL	1.156 mL
50 mM	0.0231 mL	0.1156 mL	0.2312 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481