


HPK1-IN-15

Chemical Properties

CAS No. : 2201098-03-3
 Formula: C₂₄H₂₁ClF₃N₅
 Molecular Weight: 471.91
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	HPK1-IN-15 is a selective and effective inhibitor of HPK1. Hematopoietic progenitor kinase 1 (HPK1), which was originally cloned from hematopoietic progenitor cells, belongs to the MAP4Ks family. HPK1-IN-15 can be used to study, prevent, or ameliorate diseases or conditions associated with HPK1 activity, such as cancer.
Targets(IC50)	Others,MAPK

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.119 mL	10.5952 mL	21.1905 mL
5 mM	0.4238 mL	2.119 mL	4.2381 mL
10 mM	0.2119 mL	1.0595 mL	2.119 mL
50 mM	0.0424 mL	0.2119 mL	0.4238 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481