

HI-346

Chemical Properties


CAS No. : 149486-68-0

Formula: C₁₄H₁₈BrN₃S

Molecular Weight: 340.28

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	HI-346 is a potent inhibitor of multidrug-resistant human immunodeficiency virus-1 (HIV-1).
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9388 mL	14.6938 mL	29.3876 mL
5 mM	0.5878 mL	2.9388 mL	5.8775 mL
10 mM	0.2939 mL	1.4694 mL	2.9388 mL
50 mM	0.0588 mL	0.2939 mL	0.5878 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481