

Compound N014-0004

## Chemical Properties


CAS No. :

Formula: C<sub>27</sub>H<sub>34</sub>O<sub>7</sub>

Molecular Weight: 470.562

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Compound N014-0004 is a useful organic compound for research related to life sciences and the catalog number is T131357.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1251 mL	10.6256 mL	21.2513 mL
5 mM	0.425 mL	2.1251 mL	4.2503 mL
10 mM	0.2125 mL	1.0626 mL	2.1251 mL
50 mM	0.0425 mL	0.2125 mL	0.425 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481