


4-Fluorophenethyl alcohol

Chemical Properties

CAS No. :	7589-27-7
Formula:	C ₈ H ₉ FO
Molecular Weight:	140.16
Storage:	Pure form: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	4-Fluorophenethyl alcohol, with CAS No. 7589-27-7, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 4-Fluorophenethyl alcohol provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.1347 mL	35.6735 mL	71.347 mL
5 mM	1.4269 mL	7.1347 mL	14.2694 mL
10 mM	0.7135 mL	3.5674 mL	7.1347 mL
50 mM	0.1427 mL	0.7135 mL	1.4269 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481