

## Compound T130623

## Chemical Properties


CAS No. :

Formula: C<sub>26</sub>H<sub>34</sub>O<sub>13</sub>

Molecular Weight: 554.55

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	2-[[1-(3,4-dihydroxy-5-methoxyphenyl)-5,7-dihydroxy-3-(hydroxymethyl)-6,8-dimethoxy-1,2,3,4-tetrahydronaphthalen-2-yl]methoxy]-4-(hydroxymethyl)oxolane-3,4-diol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8033 mL	9.0163 mL	18.0326 mL
5 mM	0.3607 mL	1.8033 mL	3.6065 mL
10 mM	0.1803 mL	0.9016 mL	1.8033 mL
50 mM	0.0361 mL	0.1803 mL	0.3607 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481