

Pomalidomide-C4-COOH

Chemical Properties

CAS No. : 2225940-48-5


Formula: C₁₈H₁₉N₃O₆

Molecular Weight: 373.36

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Pomalidomide-C4-COOH is a synthesized E3 ligase ligand-linker conjugate that incorporates the Pomalidomide based CRBN ligand and linker used in PROTAC technology.
Targets(IC50)	Others,E3 Ligase Ligand-Linker Conjugates

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6784 mL	13.3919 mL	26.7838 mL
5 mM	0.5357 mL	2.6784 mL	5.3568 mL
10 mM	0.2678 mL	1.3392 mL	2.6784 mL
50 mM	0.0536 mL	0.2678 mL	0.5357 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481