

Z30529104

## Chemical Properties


CAS No. : 2215884-34-5

Formula: C<sub>24</sub>H<sub>28</sub>N<sub>4</sub>O<sub>6</sub>S

Molecular Weight: 500.57

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Z30529104 is a Rpn11 inhibitor.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9977 mL	9.9886 mL	19.9772 mL
5 mM	0.3995 mL	1.9977 mL	3.9954 mL
10 mM	0.1998 mL	0.9989 mL	1.9977 mL
50 mM	0.040 mL	0.1998 mL	0.3995 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481