


2-FLUORO-PHENOXY)-PROPIONIC ACID

Chemical Properties

CAS No. :	17088-71-0
Formula:	C ₉ H ₉ FO ₃
Molecular Weight:	184.166
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-FLUORO-PHENOXY)-PROPIONIC ACID, with CAS No. 17088-71-0, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 2-FLUORO-PHENOXY)-PROPIONIC ACID provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	--

Solubility Information

Solubility	DMSO: 36.83 mg/mL (199.98 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.4298 mL	27.1488 mL	54.2977 mL
5 mM	1.086 mL	5.4298 mL	10.8595 mL
10 mM	0.543 mL	2.7149 mL	5.4298 mL
50 mM	0.1086 mL	0.543 mL	1.086 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481