


LP-360924

## Chemical Properties

CAS No. : 1984787-69-0  
 Formula: C<sub>26</sub>H<sub>30</sub>N<sub>10</sub>O<sub>3</sub>S  
 Molecular Weight: 562.65  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	LP-360924 is a bioactive chemical.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7773 mL	8.8865 mL	17.773 mL
5 mM	0.3555 mL	1.7773 mL	3.5546 mL
10 mM	0.1777 mL	0.8887 mL	1.7773 mL
50 mM	0.0355 mL	0.1777 mL	0.3555 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481