


3'-Deoxy-4-O-methylsappanol

Chemical Properties

CAS No. :	112408-68-1
Formula:	C17H18O5
Molecular Weight:	302.32
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	'3'-Deoxy-4-O-methylsappanol exhibits moderate cytotoxicity against cultured human melanoma HMV-II cells, with an IC50 value of 872.0 $\mu\text{M} \pm [2.4]$.
Targets(IC50)	Tyrosinase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3078 mL	16.5388 mL	33.0775 mL
5 mM	0.6616 mL	3.3078 mL	6.6155 mL
10 mM	0.3308 mL	1.6539 mL	3.3078 mL
50 mM	0.0662 mL	0.3308 mL	0.6616 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481