

3',3'''-Biapigenin

Chemical Properties


CAS No. : 151455-26-4

Formula: C30H18O10

Molecular Weight: 538.46

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	3',3'''-Biapigenin is a natural product for research related to life sciences. The catalog number is TN5500 and the CAS number is 151455-26-4.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8571 mL	9.2857 mL	18.5715 mL
5 mM	0.3714 mL	1.8571 mL	3.7143 mL
10 mM	0.1857 mL	0.9286 mL	1.8571 mL
50 mM	0.0371 mL	0.1857 mL	0.3714 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481