

Azido-PEG4-hydrazide-Boc

Chemical Properties

CAS No. : 1919045-01-4


Formula: C16H31N5O7

Molecular Weight: 405.45

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Azido-PEG4-hydrazide-Boc is a polyethylene glycol (PEG)-based bifunctional linker compound used in the synthesis of Proteolysis Targeting Chimeras (PROTACs)[1].
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two distinct ligands joined by a linker: one ligand targets an E3 ubiquitin ligase, while the other targets the desired protein. By leveraging the intracellular ubiquitin-proteasome system, PROTACs enable the selective degradation of target proteins [1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4664 mL	12.332 mL	24.664 mL
5 mM	0.4933 mL	2.4664 mL	4.9328 mL
10 mM	0.2466 mL	1.2332 mL	2.4664 mL
50 mM	0.0493 mL	0.2466 mL	0.4933 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481