


2,7-Dihydrohomoerysotrine

Chemical Properties

CAS No. : 51095-85-3
 Formula: C₂₀H₂₇NO₃
 Molecular Weight: 329.43
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	2,7-Dihydrohomoerysotrine possesses cardiodepressant activities on the heart of <i>B. glabrata</i> .
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0355 mL	15.1777 mL	30.3555 mL
5 mM	0.6071 mL	3.0355 mL	6.0711 mL
10 mM	0.3036 mL	1.5178 mL	3.0355 mL
50 mM	0.0607 mL	0.3036 mL	0.6071 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481