

Compound T129037

Chemical Properties


CAS No. :

Formula: C₃₀H₄₈O₄

Molecular Weight: 472.71

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	1-[(4E)-2,6-dihydroxy-6-methylhept-4-en-2-yl]-7-hydroxy-3a,6,6,9b,11a-pentamethyl-1H,2H,3H,3aH,3bH,4H,6H,7H,8H,9H,9aH,9bH,10H,11H,11aH-cyclopenta[a]phenanthren-10-one is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1155 mL	10.5773 mL	21.1546 mL
5 mM	0.4231 mL	2.1155 mL	4.2309 mL
10 mM	0.2115 mL	1.0577 mL	2.1155 mL
50 mM	0.0423 mL	0.2115 mL	0.4231 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481