


20,24-Epoxy-24-methoxy-23(24-25)abeo-dammaran-3-one

Chemical Properties

CAS No. :	1020074-97-8
Formula:	C ₃₁ H ₅₂ O ₃
Molecular Weight:	472.754
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	20,24-Epoxy-24-methoxy-23(24-25)abeo-dammaran-3-one is used as a reference standard compound, which can be used in relevant studies in the field of life sciences.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1153 mL	10.5764 mL	21.1528 mL
5 mM	0.4231 mL	2.1153 mL	4.2306 mL
10 mM	0.2115 mL	1.0576 mL	2.1153 mL
50 mM	0.0423 mL	0.2115 mL	0.4231 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481