

IND8

## Chemical Properties


CAS No. : 1778682-59-9

Formula: C17H19N5

Molecular Weight: 293.37

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	IND8 is a selective and potent $\alpha 7$ -nAChR agonist.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4087 mL	17.0433 mL	34.0866 mL
5 mM	0.6817 mL	3.4087 mL	6.8173 mL
10 mM	0.3409 mL	1.7043 mL	3.4087 mL
50 mM	0.0682 mL	0.3409 mL	0.6817 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481