


4-Benzylpiperidine

Chemical Properties

CAS No. :	31252-42-3
Formula:	C12H17N
Molecular Weight:	175.2701
Storage:	Pure form: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	4-Benzylpiperidine ,with CAS No. 31252-42-3, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 4-Benzylpiperidine provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.7055 mL	28.5274 mL	57.0548 mL
5 mM	1.1411 mL	5.7055 mL	11.411 mL
10 mM	0.5705 mL	2.8527 mL	5.7055 mL
50 mM	0.1141 mL	0.5705 mL	1.1411 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481