

Thalibrunine and Isothalibrunine, mixture of

Chemical Properties


CAS No. :

Formula: C₃₉H₄₄N₂O₈

Molecular Weight: 668.787

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Thalibrunine and Isothalibrunine, mixture of is a useful organic compound for research related to life sciences and the catalog number is T131583.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4952 mL	7.4762 mL	14.9524 mL
5 mM	0.299 mL	1.4952 mL	2.9905 mL
10 mM	0.1495 mL	0.7476 mL	1.4952 mL
50 mM	0.0299 mL	0.1495 mL	0.299 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481