


4',6-Dihydroxy-3'-methoxyaurone

Chemical Properties

CAS No. :	32396-80-8
Formula:	C ₁₆ H ₁₂ O ₅
Molecular Weight:	284.267
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	4',6-Dihydroxy-3'-methoxyaurone is a useful organic compound for research related to life sciences. The catalog number is T125042 and the CAS number is 32396-80-8.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.5178 mL	17.5889 mL	35.1778 mL
5 mM	0.7036 mL	3.5178 mL	7.0356 mL
10 mM	0.3518 mL	1.7589 mL	3.5178 mL
50 mM	0.0704 mL	0.3518 mL	0.7036 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481