


## 6-Chloropyridazine-3-carboxylic acid methyl ester

## Chemical Properties

CAS No. :	65202-50-8
Formula:	C <sub>6</sub> H <sub>5</sub> ClN <sub>2</sub> O <sub>2</sub>
Molecular Weight:	172.57
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Heterocyclic Compounds - Other heterocycles; Intermediates and Building Blocks - Electrophile
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.7947 mL	28.9737 mL	57.9475 mL
5 mM	1.1589 mL	5.7947 mL	11.5895 mL
10 mM	0.5795 mL	2.8974 mL	5.7947 mL
50 mM	0.1159 mL	0.5795 mL	1.1589 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481