


Deterenol HBr

Chemical Properties

CAS No. : 854179-30-9
 Formula: C₁₁H₁₈BrNO₂
 Molecular Weight: 276.17
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	Deterenol HBr is a beta-adrenoceptor agonist.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.621 mL	18.1048 mL	36.2096 mL
5 mM	0.7242 mL	3.621 mL	7.2419 mL
10 mM	0.3621 mL	1.8105 mL	3.621 mL
50 mM	0.0724 mL	0.3621 mL	0.7242 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481