


## Corynecin III

## Chemical Properties

CAS No. :	18048-95-8
Formula:	C <sub>13</sub> H <sub>18</sub> N <sub>2</sub> O <sub>5</sub>
Molecular Weight:	282.29
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Corynecin III is a chloramphenicol-like antibiotic originally isolated from Corynebacterium. It inhibits the growth of Gram-positive and Gram-negative bacteria, including <i>S. faecalis</i> , <i>S. aureus</i> , <i>B. subtilis</i> , <i>E. coli</i> , <i>P. aeruginosa</i> , <i>P. vulgaris</i> , <i>S. sonnei</i> , <i>S. typhosa</i> , and <i>K. pneumoniae</i> (MICs = 2.6-83 µg/ml).
Targets(IC50)	Others,Antibacterial,Antibiotic

## Solubility Information

Solubility	DMSO: Soluble (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.5425 mL	17.7123 mL	35.4246 mL
5 mM	0.7085 mL	3.5425 mL	7.0849 mL
10 mM	0.3542 mL	1.7712 mL	3.5425 mL
50 mM	0.0708 mL	0.3542 mL	0.7085 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481